3Shape E3

LAB SCANNER | MAINTENANCE SCHEDULE

The3ShapeE3sitsmid-rangeintheESeries,balancing speed and accuracy for general restorative and implantwork.Maintenancemirrorstherestoftheseries, with calibration discipline being the key to long-term accuracy.

DAILY

- Wipe the scanner exterior and chamber with a soft, lint-free cloth before the day's first scan.
- Inspect models for moisture, dust, or stone debris before scanning to protect internal optics.
- Confirm ScanIt Dental software connects properly before beginning a scanning session.

WEEKLY

- Clean the scan window/lens with 3Shape-approved materials, checking for smudges or dust film.
- Check the model holder/turntable mechanism for smooth, consistent movement.
- Review recent scan output for noise, gaps, or stitching errors.

MONTHLY

- Perform scanner calibration per manufacturer guidelines.
- Check for and install ScanIt Dental software updates.
- Inspect all cables and connectors for wear.

AS NEEDED

- Recalibrate immediately if scan accuracy or stitching quality degrades.
- Replace frayed or damaged cables right away.
- Contact 3Shape support or LabCare if cleaning does not resolve persistent scan artifacts.



NOWAK DIGITAL
855-669-2533