SprintRay Midas

CHAIRSIDE 3D PRINTER | **MAINTENANCE SCHEDULE**

The SprintRay Midas prints from sealed, single-use capsules using Digital Press Stereolithography, so there is no resin tank or vat to maintain. Routine care centers on the print glass and capsules.

DAILY

- Scan every capsule before placing it so the printer can verify the material, shade and expiration date.
- Make sure the cartridge template sits level on the print glass before starting a print.
- Keep hands clear of the print arm while printing. Discard each capsule after use.

WEEKLY

- Remove the template and inspect the print glass under bright light. Clean with 91% IPA on a lint-free towel, then glass cleaner.
- Let the glass dry completely before printing.
- Clean the touchscreen using Screen Cleaning Mode (Settings, Device).

MONTHLY

- Check for and install Midas software updates.
- Rotate capsule stock and check expiration dates. Keep capsules sealed until use.
- Clean and inspect your post-processing equipment per its instructions.

AS NEEDED

- Clean the print glass any time you see smudges, resin or debris.
- Never spray liquids directly onto the printer. Apply cleaners to a towel first.
- If restorations fail or fit changes after cleaning, call Nowak support.



NOWAK DIGITAL
855-669-2533